



## ***DNA Extraction***

The purpose of this protocol is to gather a number (in the millions!) of cells, break open the cells and nuclei, release the DNA, and cause it to stick together so you can see it!

### **Materials**

Lysis buffer  
Sterile 15 mL tube  
Isopropanol

### **Protocol:**

#### **COLLECT AND BREAK OPEN CELLS**

1. Each person in the group labels an empty sterile 15 ml tube with their name.  
Go to the water fountain located outside the laboratory and add 3 mL of water to the tube.
2. GENTLY chew on the inside of your cheek for 30 seconds. (**Do Not** draw blood.) This removes enough cells for you from which to extract DNA.
3. Pour the 3 ml of water in your mouth and rinse for 30 seconds (like mouthwash).
4. Spit water back into tube. Return to the laboratory.
5. Add 2 ml of Lysis buffer to your 15 ml tube containing “rinse” water.
6. Put in 50°C waterbath for 10 minutes.

#### **PRECIPITATION OF DNA**

1. Gently add 10 ml of ice cold isopropanol to your 15ml tube.
2. Observe 2 layers (water and isopropanol).
3. Leave in rack for 5 minutes while DNA precipitates.
4. Turn over 5 times to mix contents.
5. Stringy white/clear material is your DNA!